

Capsule Deduster VCD 100



Compact

High elevation capacity

Integrated metal detector

Empty Capsule Elimination


Efficient

Easy cleaning solution



Capsule Deduster VCD 100

Rep.	Désignation	Matière
01	Dedusting screw	Polyamide
02	Empty Capsule Eliminator ECE	St. Steel 316L
03	Metal detector CEIA	St. Steel 316L
04	Pneumatic Control Panel	St. Steel
05	Control Panel CEIA	St. Steel
06	Capsule inlet and Dedusting tube	St. Steel 316L
07a	Control Panel Danfoss	
07b	HMI Magelis SCHNEIDER	
08	Gearmotor SOTIC	
09	Speed variator Danfoss	
10	Quick dismantling triclamp PTFE	St. Steel 316L



Metal Detector CEIA



Empty Capsule Eliminator ECE



Level scanning System LRD



Feeding hopper TVD



Capsule Diverter RGV



Receiving hopper TVR



TECHNIlab®

30 bis, rue du bois chaland, CE 2948 - 91029 EVRY Cedex France

Tél : + 33 (0) 1 60 86 96 48 - Fax : + 33 (0) 1 60 86 90 01

www.technilab.fr